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1. HPLC Columns, Theory, Technology and Practice. Uwe Neue, Wiley-VCH 1997.
2. Jerkovitch AD, Mellors JS, and Jorgenson JW. LCGC North America, 2003; 21: 7
3. Wu N, Lippert JA, and Lee ML. J. Chromatogr. 2001; 911: 1
4. Wyndham KD, O'Gara JE, et al. Anal. Chem. 2003; 75: 6781-6788
5. Application Notes at [www.waters.com](http://www.waters.com) (search by Library No.)

Simultaneous Quantification of Chemically Diverse Compounds Using ESCI™ Multi-Mode Ionization on a Waters ACQUITY UPLC™-Quattro Premier™

Author: Kate Yu1, Peter Alden1, Li Di2, Susan Q. Li2, Ed Kerns2, Source: Waters Applications Note / Year: 2005 / Volume: / Page: 4 pp. Library No.: 720001097EN

HT Quantitative Analysis for A Drug Mixture by LC/MS/MS: UPLC™/MS/MS and HPLC/MS/MS Compared. Author: Kate Yu, David Little, Rob Plumb, Source: Waters Applications Note / Year: 2005 / Volume: February / Page: 4 pp. Library No.: 720001120EN

Linearity of an ACQUITY UPLC™/Quattro Premier™ System for the Analysis of 17-Hydroxyprogesterone in Protein-Precipitated Plasma, Author: Tabisam Khan and Lisa J Calton, Source: Waters Applications Note / Year: 2005 / Volume: / Page: 3 pp. Library No.: 720001135EN

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
based Bioanalytical Assays Using UPLC™, Author: Nicholas Ellor, Frances Gorycki, Chung-Ping Yu, Source: Waters Applications Note / Year: 2005 / Volume: February / Page: 2 pp. Library No.: 720001122EN

A New Paradigm for Metabolism Studies: UPLC/Q-ToF, Author: Jose Castro-Perez, Robert Plumb, Jennifer Granger, Waters Corporation, Source: Waters Applications Note / Year: 2004 / Volume: / Page: 6pp. Library No.: 720000953EN

UPLC-*oa*TOF MS with One-Minute Separation Times Applied to a Metabonomics 90-Day Tox Study, Author: Jennifer H. Granger and Robert S. Plumb, Waters Corporation; John N. Haselden, Maria L. Beaumont and Mark Hodson, GlaxoSmithKline, Source: Waters Applications Note / Year: 2004 / Volume: / Page: 5 pp. Library No.: 720001054EN

Non-Mammalian Metabonomics: An Analysis of Coffee Bean Extracts, Author: Jennifer H. Granger, Robert S. Plumb, and James N. Willis, Waters Corporation, Source: Waters Applications Note / Year: 2004 / Volume: / Page: 3 pp. Library No.: 720001058EN

Metabonomics Analysis of Zucker Rat Urine using UPLC/MS(TOF): A Time-Course Study, Author: Prof. Ian D. Wilson, Rebecca Williams, AstraZeneca; Jennifer H. Granger and Robert S. Plumb, Waters Corporation, Source: Waters Applications Note / Year: 2004 / Volume: / Page: 4 pp. Library No.: 720001057EN

Application of Ultra Performance Liquid Chromatography for LC/MS-Based Metabonomics, Author: Prof. Ian Wilson, AstraZeneca; Jennifer Granger, Robert Plumb, Waters, Source: / Year: 2004 / Volume: / Page: Library No.: 720000866EN


The Metabolism of Acetaminophen: Harnessing the Power of UPLC/MS: Jennifer H. Granger and Robert S. Plumb, Waters Corporation, Source: Waters Applications Note / Year: 2004 / Volume: / Page: 2 pp. Library No.: 720001059EN

ACQUITY UPLC for the Rapid Analysis of Soft Drinks, Author: Andrew Aubin, Waters Corporation, Year: 2004 / Volume: / Page: 4 pp. Library No.: 720001053EN

Maximizing Chromatographic Resolution of Metabolites Using UPLC, Author: Jose Castro-Perez, Robert Plumb, Jennifer Granger, Source: LC/GC North America Applications Notebook / Year: 2004 / Volume: September / Page: 1 Library No.: 720001016EN

Comparison of the Resolution and Sensitivity of UPLC and Monolithic Columns for Metabolite Profiling, Author: Robert S. Plumb, Jose Castro-Perez, Kelly A. Johnson, and Michael D. Jones, Waters Corporation, Source: Waters Applications Note / Year: 2004 / Volume: / Page: 2 pp. Library No.: 720001064EN

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
HPLC conditions:

HPLC System: Waters® Alliance® HT 2795

Column: Waters Symmetry® C18 (250 mm x 4.6 mm, 5 mm particle size) with guard column (20 mm x 3.9 mm, 5 mm particle size).

Mobile phase: MeCN (B) : H<sub>2</sub>O (0.2% HCOOH) (A) • Gradient: 0-10 min: 15% B; 10-40 min: 15-30% B; 40-50 min: 30-15% B

Column temperature: 35°C; Flow: 1 mL/min - split 1:4

UPLC conditions:

System: Waters ACQUITY UPLC™

Column: Waters ACQUITY UPLC™ BEH C18 (100 mm x 2.1 mm, 1.7 mm particle size)

Mobile phase: MeCN (B) : H<sub>2</sub>O (0.2% HCOOH) (A)

Gradient: 0-1 min: 5% B; 1-10 min: 20% B; 10-13 min: 20% B; 13-14 min: 30% B; 14-19 min: 30% B; 19.1-27 min 5% B

Column temperature: 40°C; Flow: 0.6 mL/min

Mass Spectrometry conditions:

MS= LCT Premier™ *oa*-TOF

Capillary Voltage: 3000 V (+VE) and 2600 V (-VE)

Ionisation mode: Positive and Negative electrospray

Resolution: 5500 FWHM (V mode) and 12000 FWHM (W mode)

Reference Lock mass: Leucine Enkephalin [M+H]<sup>+</sup>= 556.2771

Leucine Enkephalin [M-H]<sup>-</sup>= 554.2615

LockSpray switch time: 10 spectra HPLC

30 spectra ACQUITY UPLC™

Acquisition time: 1 spectra/second HPLC

1 spectra/ 0.2 seconds ACQUITY UPLC™